

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rabbit anti Bacterial Pyruvate oxidase, Clone: [Polyclonal], Monoclonal NMB-NE186/7S

|                            |                                                                                                                                                                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Rabbit anti Bacterial Pyruvate oxidase, Clone: [Polyclonal], Monoclonal                                                                                                                                                                                         |
| Biozol Catalog Number      | NMB-NE186/7S                                                                                                                                                                                                                                                    |
| Supplier Catalog Number    | NE186/7S                                                                                                                                                                                                                                                        |
| Alternative Catalog Number | NMB-NE186/7S                                                                                                                                                                                                                                                    |
| Manufacturer               | NordicMubio                                                                                                                                                                                                                                                     |
| Host                       | Rabbit                                                                                                                                                                                                                                                          |
| Category                   | Antikörper                                                                                                                                                                                                                                                      |
| Product Description        | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...   |
| Clonality                  | Monoclonal                                                                                                                                                                                                                                                      |
| Clone Designation          | [Polyclonal]                                                                                                                                                                                                                                                    |
| Buffer                     | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstituted by adding 1.0 ml sterile distilled |
| Source                     | Pyruvate oxidase isolated and purified from bacteria. Freund's complete adjuvant is used in the first step of the immunization procedure.                                                                                                                       |

|                   |                                                                                                                                                                                                          |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formula           | Purified hyperimmune rabbit IgG lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | Indirect immunofluorescence,ELISA,Dot blot,Immunoblotting.                                                                                                                                               |