

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Rabbit anti Saccharomyces 6-Phosphogluconic dehydrogenase, Clone: [Polyclonal], Monoclonal NMB-NE128/PAB**

|                            |                                                                                                                                                                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Rabbit anti Saccharomyces 6-Phosphogluconic dehydrogenase, Clone: [Polyclonal], Monoclonal                                                                                                                                                                      |
| Biozol Catalog Number      | NMB-NE128/PAB                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | NE128/PAb                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | NMB-NE128/PAB                                                                                                                                                                                                                                                   |
| Manufacturer               | NordicMubio                                                                                                                                                                                                                                                     |
| Host                       | Rabbit                                                                                                                                                                                                                                                          |
| Category                   | Antikörper                                                                                                                                                                                                                                                      |
| Product Description        | The reagents were evaluated for potency, purity and specificity using most or all of the following techniques: immunoelectrophoresis, cross-immunoelectrophoresis, single radial immunodiffusion (Ouchterlony), block titration, ELISA, immunoblotting and...   |
| Clonality                  | Monoclonal                                                                                                                                                                                                                                                      |
| Clone Designation          | [Polyclonal]                                                                                                                                                                                                                                                    |
| Buffer                     | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. It is reconstitute |
| Source                     | 6-Phosphogluconic dehydrogenase isolated and purified from yeast. Freund's complete adjuvant is used in the first step of the immunization procedure.                                                                                                           |

|                   |                                                                                                                                                                                                                                              |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formula           | Purified hyperimmune rabbit IgG antibodies lyophilised from a solution in phosphate buffered saline (PBS, pH 7.2), stabilized with dextran. No preservative added, as it may interfere with the antibody activity. No foreign protein added. |
| Application Notes | Immunoprecipitation.                                                                                                                                                                                                                         |