

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Llama IgG Affinity Resin (For the removal of cross-reactivity to llama immunoglobulins), Unconjugated IMU-GE-110

Article Name	Llama IgG Affinity Resin (For the removal of cross-reactivity to llama immunoglobulins), Unconjugated
Biozol Catalog Number	IMU-GE-110
Supplier Catalog Number	GE-110
Alternative Catalog Number	IMU-GE-110
Manufacturer	ImmunoReagents
Category	Antikörper
Conjugation	Unconjugated
Concentration	8-10 mg protein/ml resin
Buffer	10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Form	Liquid solution containing 4% agarose beads