

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Opal Prestained Protein Standard 3.5-245kDa ROC-MB-211-0500

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Opal Prestained Protein Standard 3.5-245kDa                                                                                                                                                                                                                   |
| Biozol Catalog Number      | ROC-MB-211-0500                                                                                                                                                                                                                                               |
| Supplier Catalog Number    | MB-211-0500                                                                                                                                                                                                                                                   |
| Alternative Catalog Number | ROC-MB-211-0500                                                                                                                                                                                                                                               |
| Manufacturer               | Rockland Immunochemicals                                                                                                                                                                                                                                      |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                              |
| Application                | SDS-PAGE, WB                                                                                                                                                                                                                                                  |
| Conjugation                | Unconjugated                                                                                                                                                                                                                                                  |
| Concentration              | Neat                                                                                                                                                                                                                                                          |
| Buffer                     | See application note.                                                                                                                                                                                                                                         |
| Form                       | Liquid (in 1x Loading Buffer)                                                                                                                                                                                                                                 |
| Application Notes          | Opal Prestained Protein Standard 3.5-245kDa is supplied in gel loading buffer and is ready to use. Do not heat, dilute, add reducing agent before loading. Suggested load volumes: 3 µl or 5 µl per loading for clear visualization during electrophoresis on |