

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Sapphire Prestained Protein Standard 10-180kDa ROC-MB-212-0500

Article Name	Sapphire Prestained Protein Standard 10-180kDa
Biozol Catalog Number	ROC-MB-212-0500
Supplier Catalog Number	MB-212-0500
Alternative Catalog Number	ROC-MB-212-0500
Manufacturer	Rockland Immunochemicals
Category	Proteine/Peptide
Application	SDS-PAGE, WB
Conjugation	Unconjugated
Buffer	See application note.
Form	Liquid (in 1x Loading Buffer)
Application Notes	Sapphire Prestained Protein Standard 10-180kDa is supplied in gel loading buffer and is ready to use. Do not heat, dilute, add reducing agent before loading. Suggested load volumes: 3 µl or 5 µl per loading for clear visualization during electrophoresis