

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Topoisomerase II DNA Decatenation Assay Kit Plus-100 (enzyme included) PFD-HDC100KE

|                            |                                                                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human Topoisomerase II DNA Decatenation Assay Kit Plus-100 (enzyme included)                                                                                                  |
| Biozol Catalog Number      | PFD-HDC100KE                                                                                                                                                                  |
| Supplier Catalog Number    | HDC100KE                                                                                                                                                                      |
| Alternative Catalog Number | PFD-HDC100KE                                                                                                                                                                  |
| Manufacturer               | ProFoldin                                                                                                                                                                     |
| Category                   | Molekularbiologie                                                                                                                                                             |
| Product Description        | This product includes the reaction buffer, concatenated DNA, human topoisomerase II alpha, spin columns and loading solution for 100 assays of DNA decatenation reactions.... |