

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Label IT Nucleic Acid Labeling Kit, Fluorescein MIR-3200

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Label IT Nucleic Acid Labeling Kit, Fluorescein                                                                                                                                                                                                               |
| Biozol Catalog Number      | MIR-3200                                                                                                                                                                                                                                                      |
| Supplier Catalog Number    | MIR-3200                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | MIR-3200,MIR3200,MIR 3200                                                                                                                                                                                                                                     |
| Manufacturer               | Mirus Bio                                                                                                                                                                                                                                                     |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Product Description        | With Label IT fluorophores and labels, you can quickly and easily label any purified single or double-stranded nucleic acid thats at least 5 nucleotides long. Simply choose the right Label IT reagent, add it to the nucleic acid, incubate for one hour... |