

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Monosan HIER Tris-EDTA Buffer pH 8,0(10x), for 5.000 ml MON-APP162**

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Monosan HIER Tris-EDTA Buffer pH 8,0(10x), for 5.000 ml                                                                                                                                                                                                       |
| Biozol Catalog Number      | MON-APP162                                                                                                                                                                                                                                                    |
| Supplier Catalog Number    | MON-APP162                                                                                                                                                                                                                                                    |
| Alternative Catalog Number | MON-APP162                                                                                                                                                                                                                                                    |
| Manufacturer               | Monosan                                                                                                                                                                                                                                                       |
| Category                   | Sonstiges                                                                                                                                                                                                                                                     |
| Application                | IF, IHC-Fr, IHC-P                                                                                                                                                                                                                                             |
| Product Description        | Immunohistochemical staining procedures consist of sequential incubation steps with blocking solutions, antibodies and secondary reagents, enzymes and chromogenic substrates carried out on tissue sections. These tissue sections are mostly prepared ou... |