

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Glucose Transporter 3 (GLUT-3 / GLUT3) ELISA KBS-KBH1434

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Glucose Transporter 3 (GLUT-3 / GLUT3) ELISA                                                                                                                                                                                                    |
| Biozol Catalog Number      | KBS-KBH1434                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | KBH1434                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | KBS-KBH1434-1                                                                                                                                                                                                                                                 |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Human                                                                                                                                                                                                                                                         |
| Product Description        | The GENLISA Human Glucose Transporter 3 (GLUT-3 / GLUT3) ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistency - Accuracy and Precision Validat... |
| Range                      | 0.78 - 50 ng/ml                                                                                                                                                                                                                                               |
| Sensitivity                | 0.31 ng/ml                                                                                                                                                                                                                                                    |
| Samples                    | Serum, Plasma, Cell Culture Supernatant and other biological samples. Validated for citrated/EDTA Plasma only.                                                                                                                                                |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |