

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Hot Start E-Taq DNA Polymerase, Glycerol-Free (5 U/ μ L) SCX-TQ08-E21B-500

| | |
|----------------------------|---|
| Article Name | Hot Start E-Taq DNA Polymerase, Glycerol-Free (5 U/ μ L) |
| Biozol Catalog Number | SCX-TQ08-E21B-500 |
| Supplier Catalog Number | TQ08-E21B-500 |
| Alternative Catalog Number | SCX-TQ08-E21B-500 |
| Manufacturer | SignalChem Diagnostics |
| Category | Molekularbiologie |
| Product Description | Hot Start E-Taq DNA Polymerase, Glycerol-Free (5 U/ μ L) features dual antibody-based blocking with inhibited 5a+3 polymerase and 5a+3 exonuclease activities.... |