

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Cluster of Differentiation 44 (CD44) ELISA KBS-KB3071

|                            |                                                                                                                                                                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Cluster of Differentiation 44 (CD44) ELISA                                                                                                                                                                                                     |
| Biozol Catalog Number      | KBS-KB3071                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | KB3071                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | KBS-KB3071-1                                                                                                                                                                                                                                                 |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                          |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                  |
| Species Reactivity         | Human                                                                                                                                                                                                                                                        |
| Product Description        | The GENLISA Human CD44 ELISA includes features like: - Ready to use protocol with break-apart wells for ease of use - Standardisation and High Reproducibility - Lot to Lot Consistency - Accuracy and Precision Validated against seven points for a "GO... |
| Range                      | 125 - 8000 pg/ml                                                                                                                                                                                                                                             |
| Sensitivity                | 75 pg/ml                                                                                                                                                                                                                                                     |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                           |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                   |