

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **GENLISA Human Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2A)** **KBS-KBH10711**

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2A)                                                                                                                                                                                   |
| Biozol Catalog Number      | KBS-KBH10711                                                                                                                                                                                                                                                  |
| Supplier Catalog Number    | KBH10711                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | KBS-KBH10711-1                                                                                                                                                                                                                                                |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                              |
| Species Reactivity         | Human                                                                                                                                                                                                                                                         |
| Product Description        | The GENLISA Human Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2A) includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precisio... |
| Range                      | 0 - 1000 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                    |
| Sensitivity                | 30 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                          |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                            |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |