

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### NexaTag Canine IGFBP-2 qIPCR ELISA Kit RAY-BIQC-IGFBP2-2

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | NexaTag Canine IGFBP-2 qIPCR ELISA Kit                                                                                                                                                                                                                        |
| Biozol Catalog Number      | RAY-BIQC-IGFBP2-2                                                                                                                                                                                                                                             |
| Supplier Catalog Number    | BIQC-IGFBP2-2                                                                                                                                                                                                                                                 |
| Alternative Catalog Number | RAY-BIQC-IGFBP2-2                                                                                                                                                                                                                                             |
| Manufacturer               | RayBiotech                                                                                                                                                                                                                                                    |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Canine                                                                                                                                                                                                                                                        |
| Product Description        | The bead-based NexaTag Canine IGFBP-2 qIPCR ELISA Kit merges the specificity of sandwich-based ELISA, the sensitivity of real-time PCR, and the easy manipulation of microspheres to detect even the lowest levels of IGFBP-2 in serum, plasma, and cell c... |
| Form                       | 96-well Microplate                                                                                                                                                                                                                                            |
| Samples                    | Cell Culture Supernatants, Plasma, Serum                                                                                                                                                                                                                      |
| Target                     | IGFBP-2 (IGFBP2)                                                                                                                                                                                                                                              |
| Application Notes          | qPCR                                                                                                                                                                                                                                                          |