

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Chlamydia Trachomatis HSP70 (549-660 a.a.) Recombinant RAY-228-10232-2

Article Name	Chlamydia Trachomatis HSP70 (549-660 a.a.) Recombinant
Biozol Catalog Number	RAY-228-10232-2
Supplier Catalog Number	228-10232-2
Alternative Catalog Number	RAY-228-10232-2
Manufacturer	RayBiotech
Category	Proteine/Peptide
Species Reactivity	Bacteria
Expression System	E. coli
Purity	Chlamydia HSP70 protein was purified by proprietary chromatographic technique.
Formula	50mM Tris-HCl, pH 8.0, 8M Urea, 60mM NaCl and 50% glycerol.
Application Dilute	Chlamydia HSP70 protein is suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.