

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

MabSelect SuRe pcc protein A resin CRM-MA-331C

Article Name	MabSelect SuRe pcc protein A resin
Biozol Catalog Number	CRM-MA-331C
Supplier Catalog Number	Ma-331C
Alternative Catalog Number	CRM-MA-331C-25
Manufacturer	Creative BioMart
Category	Proteine/Peptide
Application Notes	<p>The MabSelect family, comprising MabSelect, MabSelect Xtra, MabSelect SuRe, and MabSelect SuRe LX chromatography media, was developed for high-capacity MAb capture at process scale. To enable the high mass transfer required for continuous chromatography, MabSelect SuRe pcc was further developed from MabSelect SuRe LX. The high dynamic binding capacity (DBC) for MAbs at short residence time contributes to increased productivity by reducing total time required for column load operations. Comparative studies with MabSelect SuRe LX show that higher capacities can be achieved at shorter residence times with MabSelect SuRe pcc.</p>