

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Link-Dtech 545-Tb (Green) - For LFIA, Bright-Dtech 545 PDT-LDT545-LFIA

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Link-Dtech 545-Tb (Green) - For LFIA, Bright-Dtech 545                                                                                                                                                                                                        |
| Biozol Catalog Number      | PDT-LDT545-LFIA                                                                                                                                                                                                                                               |
| Supplier Catalog Number    | LDT545 - LFIA                                                                                                                                                                                                                                                 |
| Alternative Catalog Number | PDT-LDT545-LFIA-1                                                                                                                                                                                                                                             |
| Manufacturer               | Poly-Dtech                                                                                                                                                                                                                                                    |
| Category                   | Sonstiges                                                                                                                                                                                                                                                     |
| Application                | LF                                                                                                                                                                                                                                                            |
| Conjugation                | Bright-Dtech 545                                                                                                                                                                                                                                              |
| Product Description        | Link-Dtech 545-Tb (Green) kit is coupled with Bright-Dtech 545 for obtaining a labeling of antibodies, peptides, oligonucleotides and proteins. The kit is perfect for the need of small-volume fast labeling, for TRF assay development (FRET, FLISA, WB,... |
| Range                      | 1x10 <sup>-8</sup> M                                                                                                                                                                                                                                          |
| Sensitivity                | 1x10 <sup>-13</sup> M or Possibility &lt; 1 pg/mL                                                                                                                                                                                                             |
| Form                       | Liquid                                                                                                                                                                                                                                                        |
| Samples                    | Serum, Plasma, Cell supernatant, urine, other biological fluids                                                                                                                                                                                               |
| Application Notes          | 4-8 C after coupling                                                                                                                                                                                                                                          |