

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### BA (Butyric Acid) ELISA Kit EKL-EU21247

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | BA (Butyric Acid) ELISA Kit                                                                                                                                                                                                                                   |
| Biozol Catalog Number      | EKL-EU21247                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | EU21247                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-EU21247-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | All                                                                                                                                                                                                                                                           |
| Product Description        | This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with BA. During the reaction, BA in the sample or standard competes with a fixed amount of BA on the solid phase supporter for sites o... |
| Range                      | 0.124-10ng/mL                                                                                                                                                                                                                                                 |
| Sensitivity                | 0.056ng/mL                                                                                                                                                                                                                                                    |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | BA                                                                                                                                                                                                                                                            |