

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human GIPR (GFP) HEK293 Reporter Cell Preis auf Anfrage ABS-CHEK-ATF268-2VIALS

|                            |                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------|
| Article Name               | Human GIPR (GFP) HEK293 Reporter Cell Preis auf Anfrage                                      |
| Biozol Catalog Number      | ABS-CHEK-ATF268-2VIALS                                                                       |
| Supplier Catalog Number    | CHEK-ATF268-2Vials                                                                           |
| Alternative Catalog Number | ABS-CHEK-ATF268-2VIALS                                                                       |
| Manufacturer               | AcroBiosystems                                                                               |
| Host                       | Human                                                                                        |
| Category                   | Zellen/Zellkultur                                                                            |
| Species Reactivity         | Human                                                                                        |
| Product Description        | Human GIPR (GFP) HEK293 Reporter Cell is expressed from Human HEK293. - Zellen/Zellkultur... |
| Tag                        | GFP                                                                                          |
| Form                       | Liquid                                                                                       |
| Target                     | GIPR                                                                                         |
| Application Notes          | Screen for agonists that can bind and activate GIPR                                          |