

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Anti-Mouse IgG-coupled Magnetic Beads (recommended for MPCLIA), Human ABS-MPC-A003-100MG

Article Name	Anti-Mouse IgG-coupled Magnetic Beads (recommended for MPCLIA), Human
Biozol Catalog Number	ABS-MPC-A003-100MG
Supplier Catalog Number	MPC-A003-100mg
Alternative Catalog Number	ABS-MPC-A003-100MG
Manufacturer	AcroBiosystems
Host	Human
Category	Sonstiges
Species Reactivity	Rabbit
Product Description	AntiMouse IgGcoupled Magnetic Beads (recommended for MPCLIA) is expressed from Rabbit HEK293. - Sonstiges...
Form	Powder
Target	Mouse IgG
Application Notes	The Anti-Mouse IgG-coupled Magnetic Beads is used to capture the mouse IgG, it can combination with Acridinium ester markers in chemiluminescence technology, The Acridinium ester markers such as Streptavidin-Acridinium ester can capture the biotinylated proteins or molecules, this allows detection of antigen and antibody binding or antibody screening.