

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Extracellular OCR Plate Assay Kit DOJ-E297-12

|                            |                                                                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Extracellular OCR Plate Assay Kit                                                                                                                                                                   |
| Biozol Catalog Number      | DOJ-E297-12                                                                                                                                                                                         |
| Supplier Catalog Number    | E297-12                                                                                                                                                                                             |
| Alternative Catalog Number | DOJ-E297-12                                                                                                                                                                                         |
| Manufacturer               | Dojindo                                                                                                                                                                                             |
| Category                   | Kits/Assays                                                                                                                                                                                         |
| Product Description        | Applicable to regular fluorescent plate reader with temperature-controlled incubation No need for an expensive instrument, special medium, and plates All-in-One Kit with OCR calculation Sheets... |