

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Glucose Uptake Plate Assay Kit DOJ-UP08-10

|                            |                                                                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Glucose Uptake Plate Assay Kit                                                                                                                                                        |
| Biozol Catalog Number      | DOJ-UP08-10                                                                                                                                                                           |
| Supplier Catalog Number    | UP08-10                                                                                                                                                                               |
| Alternative Catalog Number | DOJ-UP08-10                                                                                                                                                                           |
| Manufacturer               | Dojindo                                                                                                                                                                               |
| Category                   | Kits/Assays                                                                                                                                                                           |
| Product Description        | 1. Glucose uptake capacity can be easily measured with a microplate reader, 2. Non-wash type that does not require any washing, 3. Simple operation: only requires adding reagents... |