

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EasyStep Human VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit EBT-ELK088ES

Article Name	EasyStep Human VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit
Biozol Catalog Number	EBT-ELK088ES
Supplier Catalog Number	ELK088ES
Alternative Catalog Number	EBT-ELK088ES-96, EBT-ELK088ES-48, EBT-ELK088ES-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Human
Concentration	10 µg/mL
Range	0.16-10 µg/mL
Sensitivity	32.1 ng/mL
Samples	Serum, plasma

Application Notes	Assay Type: Sandwich. Assay length: 1.5h. Research Area: CD & Adhesion molecule,Tumor immunity,Infection immunity,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human VE-Cadherin, and the Human VE-Cadherin standard plate wells that pre-coated using protein-related techniques are provided separately. Standard/Sample Diluent Buffer or samples are added to the appropriate microtiter plate wells ,then added a HRP-conjugated antibody specific to Human VE-Cadherin. After TMB substrate solution is added, only those wells that contain Human VE-Cadherin and HRP-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Human VE-Cadherin in the samples is then determined by comparing the OD of the samples to the standard curve
-------------------	--