

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

SuperKine(TM) Enhanced Antifade Mounting Medium with DAPI ASC-BMU107-EN

Article Name	SuperKine(TM) Enhanced Antifade Mounting Medium with DAPI
Biozol Catalog Number	ASC-BMU107-EN
Supplier Catalog Number	BMU107-EN
Alternative Catalog Number	ASC-BMU107-EN-10, ASC-BMU107-EN-50
Manufacturer	Abbkine Scientific
Category	Zellen/Zellkultur
Product Description	SuperKine(TM) Enhanced Antifade Mounting Medium with DAPI is compatible with a wide range of fluorescent dyes and protects the entire visible and infrared spectrum from fluorescence signal quenching....
Tag	SuperKine(TM) Enhanced Antifade Mounting Medium with DAPI
Application Notes	SuperKine(TM) Enhanced Antifade Mounting Medium with DAPI is used to seal fluorescent tissue and cell samples with DAPI and unique anti-fluorescence quencher ingredients and compatible with a variety of fluorescent dyes to protect the entire visible and infrared spectrum and prevent fluorescence signal quenching. The original luminescence intensity can be maintained more than 2 weeks when the sample is stored at 4°C or -20°C.