

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Vasorin (VASN) ELISA Kit ASC-KTE60072

Article Name	Human Vasorin (VASN) ELISA Kit
Biozol Catalog Number	ASC-KTE60072
Supplier Catalog Number	KTE60072
Alternative Catalog Number	ASC-KTE60072-48, ASC-KTE60072-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Vasorin (VASN) ELISA Kit employs a two-site sandwich ELISA to quantitate VASN....
Range	Please inquire
Sensitivity	Please inquire
Tag	VASN
NCBI	114990
UniProt	Q6EMK4
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Human VASN (VASN) ELISA Kit employs a two-site sandwich ELISA to quantitate VASN in samples. An antibody specific for VASN has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any VASN present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for VASN is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of VASN bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	--