

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Taste receptor type 2 member 19 (TAS2R19) ELISA Kit ASC-KTE60491

|                            |                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human Taste receptor type 2 member 19 (TAS2R19) ELISA Kit                                                                  |
| Biozol Catalog Number      | ASC-KTE60491                                                                                                               |
| Supplier Catalog Number    | KTE60491                                                                                                                   |
| Alternative Catalog Number | ASC-KTE60491-48, ASC-KTE60491-96                                                                                           |
| Manufacturer               | Abbkine Scientific                                                                                                         |
| Category                   | Kits/Assays                                                                                                                |
| Application                | ELISA                                                                                                                      |
| Species Reactivity         | Human                                                                                                                      |
| Product Description        | This Human Taste receptor type 2 member 19 (TAS2R19) ELISA Kit employs a two-site sandwich ELISA to quantitate TAS2R19.... |
| Range                      | Please inquire                                                                                                             |
| Sensitivity                | Please inquire                                                                                                             |
| Tag                        | TAS2R19                                                                                                                    |
| NCBI                       | <a href="#">259294</a>                                                                                                     |
| UniProt                    | <a href="#">P59542</a>                                                                                                     |
| Samples                    | Cell culture supernatantsSerumPlasmaOther biological fluids                                                                |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>This Human Taste receptor type 2 member 19 (TAS2R19) ELISA Kit employs a two-site sandwich ELISA to quantitate TAS2R19 in samples. An antibody specific for TAS2R19 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TAS2R19 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TAS2R19 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TAS2R19 bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|