

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (ACAP1) ELISA Kit ASC-KTE60873

Article Name	Human Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (ACAP1) ELISA Kit
Biozol Catalog Number	ASC-KTE60873
Supplier Catalog Number	KTE60873
Alternative Catalog Number	ASC-KTE60873-48, ASC-KTE60873-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (ACAP1) ELISA Kit employs a two-site sandwich ELISA to quantitate ACAP1....
Range	Please inquire
Sensitivity	Please inquire
Tag	ACAP1
NCBI	9744
UniProt	Q15027
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Human Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (ACAP1) ELISA Kit employs a two-site sandwich ELISA to quantitate ACAP1 in samples. An antibody specific for ACAP1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any ACAP1 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for ACAP1 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of ACAP1 bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	---