

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit ASC-KTE61169

Article Name	Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit
Biozol Catalog Number	ASC-KTE61169
Supplier Catalog Number	KTE61169
Alternative Catalog Number	ASC-KTE61169-48, ASC-KTE61169-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit employs a two-site sandwich ELISA to quantitate PREX2....
Range	Please inquire
Sensitivity	Please inquire
Tag	PREX2
NCBI	80243
UniProt	Q70Z35
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (PREX2) ELISA Kit employs a two-site sandwich ELISA to quantitate PREX2 in samples. An antibody specific for PREX2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any PREX2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PREX2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PREX2 bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	--