

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L) ELISA Kit ASC-KTE61487

Article Name	Human NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L) ELISA Kit
Biozol Catalog Number	ASC-KTE61487
Supplier Catalog Number	KTE61487
Alternative Catalog Number	ASC-KTE61487-48, ASC-KTE61487-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L) ELISA Kit employs a two-site sandwich ELISA to quantitate MT-ND4L....
Range	Please inquire
Sensitivity	Please inquire
Tag	MT-ND4L
NCBI	4539
UniProt	P03901
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Human NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L) ELISA Kit employs a two-site sandwich ELISA to quantitate MT-ND4L in samples. An antibody specific for MT-ND4L has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MT-ND4L present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for MT-ND4L is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MT-ND4L bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	---