

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human Beta-1,3-N-acetylglucosaminyltransferase manic fringe (MFNG) ELISA Kit ASC-KTE61635

Article Name	Human Beta-1,3-N-acetylglucosaminyltransferase manic fringe (MFNG) ELISA Kit
Biozol Catalog Number	ASC-KTE61635
Supplier Catalog Number	KTE61635
Alternative Catalog Number	ASC-KTE61635-48, ASC-KTE61635-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human Beta-1,3-N-acetylglucosaminyltransferase manic fringe (MFNG) ELISA Kit employs a two-site sandwich ELISA to quantitate MFNG....
Range	Please inquire
Sensitivity	Please inquire
Tag	MFNG
NCBI	4242
UniProt	O00587
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Human Beta-1,3-N-acetylglucosaminyltransferase manic fringe (MFNG) ELISA Kit employs a two-site sandwich ELISA to quantitate MFNG in samples. An antibody specific for MFNG has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyMFNG present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for MFNG is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MFNG bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	--