

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit ASC-KTE61728

Article Name	Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit
Biozol Catalog Number	ASC-KTE61728
Supplier Catalog Number	KTE61728
Alternative Catalog Number	ASC-KTE61728-48, ASC-KTE61728-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Human
Product Description	This Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit employs a two-site sandwich ELISA to quantitate MAPKAPK5....
Range	Please inquire
Sensitivity	Please inquire
Tag	MAPKAPK5
NCBI	8550
UniProt	Q8IW41
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Human MAP kinase-activated protein kinase 5 (MAPKAPK5) ELISA Kit employs a two-site sandwich ELISA to quantitate MAPKAPK5 in samples. An antibody specific for MAPKAPK5 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any MAPKAPK5 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for MAPKAPK5 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of MAPKAPK5 bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	---