

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Blood group antigen-binding adhesin A (Bab A) ELISA Kit ASC-KTE63021

|                            |                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human Blood group antigen-binding adhesin A (Bab A) ELISA Kit                                                                      |
| Biozol Catalog Number      | ASC-KTE63021                                                                                                                       |
| Supplier Catalog Number    | KTE63021                                                                                                                           |
| Alternative Catalog Number | ASC-KTE63021-48, ASC-KTE63021-96                                                                                                   |
| Manufacturer               | Abbkine Scientific                                                                                                                 |
| Category                   | Kits/Assays                                                                                                                        |
| Application                | ELISA                                                                                                                              |
| Species Reactivity         | Human                                                                                                                              |
| Product Description        | Human Blood group antigen-binding adhesin A (Bab A) ELISA Kit employs a two-site sandwich ELISA to quantitate Bab A in samples.... |
| Range                      | 2.5 ng/L-40 ng/L                                                                                                                   |
| Sensitivity                | 0.1 ng/L                                                                                                                           |
| Tag                        | Bab A                                                                                                                              |
| Samples                    | Cell culture supernatantsSerumPlasmaOther biological fluids                                                                        |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>Human Blood group antigen-binding adhesin A (Bab A) ELISA Kit employs a two-site sandwich ELISA to quantitate Bab A in samples. An antibody specific for Human Bab A has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any Bab A present is bound by the immobilized antibody. After removing any unbound substances, HRP-Conjugated Bab A detection antibody is added to the wells. Following a wash to remove any unbound HRP reagent, a Chromogen solution is added to the wells and color develops in proportion to the amount of Bab A bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|