

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### EliKine(TM) Human FSH ELISA Kit ASC-KTE6903

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| Article Name               | EliKine(TM) Human FSH ELISA Kit                                                               |
| Biozol Catalog Number      | ASC-KTE6903                                                                                   |
| Supplier Catalog Number    | KTE6903                                                                                       |
| Alternative Catalog Number | ASC-KTE6903-48, ASC-KTE6903-96                                                                |
| Manufacturer               | Abbkine Scientific                                                                            |
| Category                   | Kits/Assays                                                                                   |
| Application                | ELISA                                                                                         |
| Species Reactivity         | Human                                                                                         |
| Product Description        | EliKine(TM) Human FSH ELISA Kit employs a two-site sandwich ELISA to quantitate Human FSH.... |
| Range                      | 2 IU/L-70 IU/L                                                                                |
| Sensitivity                | 1.0 IU/L                                                                                      |
| Tag                        | FSHB                                                                                          |
| NCBI                       | <a href="#">2488</a>                                                                          |
| UniProt                    | <a href="#">P01225</a>                                                                        |
| Samples                    | SerumPlasmaOther biological fluids                                                            |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>EliKine(TM) Human FSH ELISA Kit employs a two-site sandwich ELISA to quantitate FSH in samples. An antibody specific for FSH has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any FSH present is bound by the immobilized antibody. After removing any unbound substances, a HRP conjugated antibody specific for FSH is added to the wells. Following a wash to remove any unbound enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of FSH bound in the initial step. The color development is stopped and the intensity of the color is measured.</p> |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|