

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Tenascin-R (TNR) ELISA Kit ASC-KTE70090

Article Name	Mouse Tenascin-R (TNR) ELISA Kit
Biozol Catalog Number	ASC-KTE70090
Supplier Catalog Number	KTE70090
Alternative Catalog Number	ASC-KTE70090-48, ASC-KTE70090-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse Tenascin-R (TNR) ELISA Kit employs a two-site sandwich ELISA to quantitate TNR....
Range	Please inquire
Sensitivity	Please inquire
Tag	TNR
NCBI	21960
UniProt	Q8BYI9
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Mouse Tenascin-R (TNR) ELISA Kit employs a two-site sandwich ELISA to quantitate TNR in samples. An antibody specific for TNR has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNR present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNR is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNR bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	--