

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Tenascin-N (TNN) ELISA Kit ASC-KTE70103

Article Name	Mouse Tenascin-N (TNN) ELISA Kit
Biozol Catalog Number	ASC-KTE70103
Supplier Catalog Number	KTE70103
Alternative Catalog Number	ASC-KTE70103-48, ASC-KTE70103-96
Manufacturer	Abbkine Scientific
Category	Kits/Assays
Application	ELISA
Species Reactivity	Mouse
Product Description	This Mouse Tenascin-N (TNN) ELISA Kit employs a two-site sandwich ELISA to quantitate TNN....
Range	Please inquire
Sensitivity	Please inquire
Tag	TNN
NCBI	329278
UniProt	Q80Z71
Samples	Cell culture supernatantsSerumPlasmaOther biological fluids

Application Notes	This Mouse Tenascin-N (TNN) ELISA Kit employs a two-site sandwich ELISA to quantitate TNN in samples. An antibody specific for TNN has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any TNN present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for TNN is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of TNN bound in the initial step. The color development is stopped and the intensity of the color is measured.
-------------------	--