
Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

**MAP Kinase p44/42 (MAPK3/MAPK1, Mitogen-activated Protein Kinase 1, MAP Kinase 1, MAPK 1, Mapk, ERT1, Extracellular Signal-regulated Kinase 2, ERK-2, MAP Kinase Isoform p42, p42-MAPK, Mitogen-activated Protein Kinase 2, MAP Kinase
USB-M2360-02G-PE**

| | |
|----------------------------|---|
| Article Name | MAP Kinase p44/42 (MAPK3/MAPK1, Mitogen-activated Protein Kinase 1, MAP Kinase 1, MAPK 1, Mapk, ERT1, Extracellular Signal-regulated Kinase 2, ERK-2, MAP Kinase Isoform p42, p42-MAPK, Mitogen-activated Protein Kinase 2, MAP Kinase |
| Biozol Catalog Number | USB-M2360-02G-PE |
| Supplier Catalog Number | M2360-02G-PE |
| Alternative Catalog Number | USB-M2360-02G-PE-100 |
| Manufacturer | US Biological |
| Host | Rabbit |
| Category | Antikörper |
| Application | FLISA, IHC, IP, WB |
| Immunogen | Synthetic peptide corresponding to aa339-353 from the C-terminal segment of rat ERK1. |
| Product Description | Mitogen activated protein kinases (MAPKs) play pivotal roles in mediating signal transduction from the cell surface to the nucleus (1) . These kinases are encoded by distinct genes and together form a family of kinases whose activation is dependent u... |
| Purity | Purified by immunoaffinity chromatography. |
| Form | Supplied as a liquid in PBS, pH 7.2. No preservatives added. Labeled with R-Phycoerythrin (PE). |