

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rat GInducible T-Cell Co Stimulator Ligand (ICOSL) ELISA Kit BYT-ORB1946908

Article Name	Rat GInducible T-Cell Co Stimulator Ligand (ICOSL) ELISA Kit
Biozol Catalog Number	BYT-ORB1946908
Supplier Catalog Number	orb1946908
Alternative Catalog Number	BYT-ORB1946908-48,BYT-ORB1946908-96
Manufacturer	Biorbyt
Category	Kits/Assays
Species Reactivity	Rat
Product Description	Rat GInducible T-Cell Co Stimulator Ligand (ICOSL) ELISA Kit...
Range	78.13-5000pg/mL
Sensitivity	46.88 pg/mL
Samples	serum, plasma, Tissue homogenate and Other biological samples
Target	ICOSLG

Application Notes

Application Notes: This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Rat ICOSL. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Rat ICOSL and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Rat ICOSL, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm. The OD value is proportional to the concentration of Rat ICOSL. You can calculate the concentration of Rat ICOSL in the samples by comparing the OD of the samples to the standard curve