

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ELISpot Pro: Human IFN-gamma (ALP) MAB-3420-2APW-10

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | ELISpot Pro: Human IFN-gamma (ALP)                                                                                                                                                                                                                            |
| Biozol Catalog Number      | MAB-3420-2APW-10                                                                                                                                                                                                                                              |
| Supplier Catalog Number    | 3420-2APW-10                                                                                                                                                                                                                                                  |
| Alternative Catalog Number | MAB-3420-2APW-10                                                                                                                                                                                                                                              |
| Manufacturer               | MABTECH                                                                                                                                                                                                                                                       |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Application                | ELISpot                                                                                                                                                                                                                                                       |
| Species Reactivity         | Human                                                                                                                                                                                                                                                         |
| Conjugation                | ALP                                                                                                                                                                                                                                                           |
| Product Description        | ELISpot Pro kit for enumeration of cells secreting IFN-gamma, The ELISpot Pro kit is Mabtechs most innovative product featuring pre-coated plates and one-step detection. The kit enables enumeration of cells secreting human IFN-gamma and contains ELIS... |
| Range                      | One-step ALP                                                                                                                                                                                                                                                  |
| NCBI                       | <a href="#">3458</a>                                                                                                                                                                                                                                          |
| Pubmed                     | <a href="https://pubmed.ncbi.nlm.nih.gov/32979941">https://pubmed.ncbi.nlm.nih.gov/32979941</a> target=_blank rel=nofollow>32979941                                                                                                                           |
| Target                     | IFN-gamma                                                                                                                                                                                                                                                     |

|                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes | <p>Complementary Products:<br/><a href="https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0">https://www.mabtech.com/products/peppool-cef-cd8-human-3616-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefta-cd4-human-3617-1-0">https://www.mabtech.com/products/peppool-cefta-cd4-human-3617-1-0</a>,<br/><a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-human-3627-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-human-3627-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0</a>,<br/><a href="https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd8-human-3637-1-0</a>, <a href="https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0">https://www.mabtech.com/products/peppool-cefras-global-cd4-and-cd8-human-3638-1-0</a>. Contents: 7-B6-1, ALP,BCIP/NBT-plus substrate for ELISpot,Pre-coated ELISpot plate, white MSIP (mAb 1-D1K),Anti-CD3 mAb (CD3-2)</p> |
|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|