

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **ViraDuctin(TM) Retrovirus Transduction Kit (200 Transductions) CBL-RV-201**

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | ViraDuctin(TM) Retrovirus Transduction Kit (200 Transductions)                                                                                                                                                                                                |
| Biozol Catalog Number      | CBL-RV-201                                                                                                                                                                                                                                                    |
| Supplier Catalog Number    | RV-201                                                                                                                                                                                                                                                        |
| Alternative Catalog Number | CBL-RV-201                                                                                                                                                                                                                                                    |
| Manufacturer               | Cell Biolabs                                                                                                                                                                                                                                                  |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Product Description        | The efficiency of retrovirus infection of target cells tends to be lower than that of adenovirus. Additives such as Polybrene can increase transduction efficiencies, but even then only a small fraction of retroviral vectors can trasduce many target c... |