

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

1 | PEB (4x) | protein extraction buffer 5 x 2 ml (4x stock) allows up to 75 isolations of plant material (using 500 µl 1x PEB for 100 mg fresh weight) or 190 isolations of algal material (using 200 µl 1x PEB for cell amounts co AGR-AS08-300

Article Name	1 PEB (4x) protein extraction buffer 5 x 2 ml (4x stock) allows up to 75 isolations of plant material (using 500 µl 1x PEB for 100 mg fresh weight) or 190 isolations of algal material (using 200 µl 1x PEB for cell amounts co
Biozol Catalog Number	AGR-AS08-300
Supplier Catalog Number	AS08-300
Alternative Catalog Number	AGR-AS08-300
Manufacturer	Agrisera
Category	Proteine/Peptide
Species Reactivity	Algae, Fungi, Plant
Product Description	PEB is an extraction buffer for disruption and solubilisation of total protein from plant tissue and algal cells. The use of the anionic detergent LDS together with the recommended procedure (combination of sonication and freeze/thaw cycles) has been...
Antibody Type	Secondary Antibody
Application Notes	protein extraction