

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit IgG anti-Horse IgG (H+L)-HRPO, MinX none DNA-SEC-182978

Article Name	Rabbit IgG anti-Horse IgG (H+L)-HRPO, MinX none
Biozol Catalog Number	DNA-SEC-182978
Supplier Catalog Number	SEC-182978
Alternative Catalog Number	DNA-SEC-182978
Manufacturer	dianova
Host	Rabbit
Category	Antikörper
Application	ELISA,IHC,WB
Species Reactivity	Equine
Immunogen	Anti-Horse IgG whole molecule was produced by repeated immunization with Horse IgG whole molecule in rabbit.
Conjugation	HRPO
Format	IgG
Target Specificity	IgG (H+L)
Cross-Adsorption (MinX)	no cross-adsorbtion
Product Description	Anti-Horse IgG Peroxidase Antibody generated in rabbit detects horse IgG. Secreted as part of the adaptive immune response by plasma B cells, immunoglobulin G constitutes 75% of serum immunoglobulins. Immunoglobulin G binds to viruses, bacteria, as w...

Clonality	Polyclonal
Concentration	2.0 mg/mL
Isotype	Ig
Buffer	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Purity	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Horse IgG coupled to agarose beads. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Rabbit Serum, Horse IgG and Horse Serum.
Form	Lyophilized
Formula	20 mM K3PO4,150 mM NaCl,pH 7,2,lyophilisate,0,01% Gentamicin
Target	Horse
Antibody Type	Secondary Antibody
Application Dilute	ELISA Dilution: 1:50,000 - 1:200,000, Immunohistochemistry Dilution: 1:500 - 1:2,500, Western Blot Dilution: 1:5,000 - 1:20,000
Application Notes	Anti-Horse IgG whole molecule is suitable for use in immunoelectrophoresis, western-blot, competitive western-blot, ELISA and competitive ELISA assays. Specific conditions for reactivity and signal detection should be optimized by the end user.