

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Chicken SAA ELISA Kit, Goat ICL-E-30SAA

|                            |                                                                              |
|----------------------------|------------------------------------------------------------------------------|
| Article Name               | Chicken SAA ELISA Kit, Goat                                                  |
| Biozol Catalog Number      | ICL-E-30SAA                                                                  |
| Supplier Catalog Number    | E-30SAA                                                                      |
| Alternative Catalog Number | ICL-E-30SAA-1                                                                |
| Manufacturer               | ICL                                                                          |
| Host                       | Goat                                                                         |
| Category                   | Kits/Assays                                                                  |
| Application                | ELISA                                                                        |
| Species Reactivity         | Gallus                                                                       |
| Product Description        | Accurate measurement of chicken serum amyloid A for inflammation studies.... |
| Range                      | 0.375 ng/ml - 24 ng/ml                                                       |
| Sensitivity                | 0.094 ng/ml                                                                  |
| Form                       | Sandwich ELISA                                                               |
| Samples                    | Plasma,Serum                                                                 |
| Target                     | Serum Amyloid A                                                              |