

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### 20-Hydroxyecdysone (20E), BLG Conjugate KER-EGA264

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | 20-Hydroxyecdysone (20E), BLG Conjugate                                                                                                                                                                                                                       |
| Biozol Catalog Number      | KER-EGA264                                                                                                                                                                                                                                                    |
| Supplier Catalog Number    | EGA264                                                                                                                                                                                                                                                        |
| Alternative Catalog Number | KER-EGA264                                                                                                                                                                                                                                                    |
| Manufacturer               | Kerafast                                                                                                                                                                                                                                                      |
| Category                   | Sonstiges                                                                                                                                                                                                                                                     |
| Product Description        | These rabbit polyclonal antibodies were generated against BSA conjugated ecdysteroids E and 20E and, in combination with BLG-conjugated ecdysteroids E and 20E, allow for a fast, accurate and low-cost means to measure titers for ecdysone (E) and 20-hy... |
| Clonality                  | Polyclonal                                                                                                                                                                                                                                                    |
| Buffer                     | Whole serum                                                                                                                                                                                                                                                   |